

#### THE GEO-MATRIX SILVERBACK DIFFERENCE





PRIME 45 TILT

FOAM 45 TILT

MAX 45 TILT



### Flexibility, Agility & Responsiveness

Blake Medical is a family-run business, and we strive to treat each client like a part of our family. If for any reason you are not happy with one of our products, please let us know and we will rectify the issue.

Our innovative engineers and committed customer service personnel are focused on ensuring a consistently positive customer experience that exceeds expectations.

### HARDWARE **IICK RELEASE**

The Quick Release mounting hardware is designed to be strong and functional. This Quick Release mounting hardware allows for infinite on chair adjustment like angle, height and depth. With a locking catch you can be assured the Silverback will never release unless you want it to.











### Product Silverback Prime Active Silverback Prime Deep Silverback Prime Extra Deep Silverback Max Active Silverback Max Deep Silverback Max Extra Deep

Blake Medical Distribution 461 Cumberland Ave | Hamilton, Ontario | L8M 2A8

Phone: 1 (866) 662-3459 Fax: 1 (905) 560-1342 Email: info@blakemedical.ca Site: www.blakemedical.ca

**Distributor Information:** 

ADP Code	ND ADP Codes		
SE0001108			
SE0001108	SEBND0015 x2		
SE0001108	SEBND0015 x3		
SE0001109			
SE0001109	SEBND0015 x2		
SE0001109	SEBND0015 x3		

'Adapting to your needs...Exceeding your expectations.'







## **Geo-Matrix Silverback Tension**

The Geo-Matrix Silverback Tension back gives you a completely custom contour. With multiple adjustable Velcro straps and quick, secure installation, this option provides an easy, flexible solution with the unique ability to accommodate other custom pieces for additional support and further customization.

The Silverback tension contains a full layer of Geo-Matrix gel for pressure reduction and highdensity foam for increased support and comfort. For added protection, the Tension also offers padding for back cane posts, and additional Geo-Matrix gel sheets are available to further increase comfort and skin protection





Standard 16"





# **Geo-Matrix Silverback Wheelchair Backs**

### Lightweight \* Spinal Support \* Maximum Comfort

Geo-Matrix Silverback wheelchair backs incorporate a vertical insert of Geo-Matrix gel to target protection to the high-risk area around the spinal column.

The spinal inserts ensure optimal comfort and skin protection from the shear that results when wheelchairs tilt and recline, reducing the risk of pressure sores from developing along the spine. Geo-Matrix gel is a mineral-oil- based, dry-polymer solid gel that stays in place

regardless of pressure and is unaffected by the pull of gravity.

unique properties of the gel allow it to move with the body, reducing shear and reactionary force.

#### **Variations Available**

The Silverback wheelchair backs are rated for up to 250 lbs. at all degrees, for years of enjoyment and use or abuse!

Where extra skin protection is necessary, the especially beneficial to address over-curvature of Silverback Max's extra layer of gel provides the spine (kyphosis). superior vibration dampening qualities, making it However, if weight is a concern, the lighter Silverback So, whether in an upright position or in full tilt, the an exceptionally comfortable and protective back Prime is also available. for power wheelchairs.

> It also delivers extra protection when the spine is subjected to added pressure from prolonged wheelchair tilt conditions.

The Silverback Max is best for most clients and is

Contact the office today to request a 30 day demo.

Need a custom back or accessory? We'll work with you to engineer a perfect solution.

Problems with a product? Live trouble-shooting & mobile pressure-mapping systems allow for onsite assessments and analysis.

## CHOOSE YOUR **ALUMINUM SHELL**

**Deep Contour** 4" Depth (100mm)



Active Contour

2.5" Depth (63.5mm)





# DETAILS

Features	Geo - Matrix Prime	Geo - Matrix Max
Three contour options for the perfect amount of support: active (2.5" depth), deep (4" depth) and extra deep (6" depth)	•	•
Three height options: standard (16"), tall (18") and extra tall (20")	٠	٠
Three width options: 16", 18" and 20"	•	•
Quick-release mounting hardware for easy installation and adjustment without having to remove the chair back or reposition the user; strong, lightweight and long-lasting	on 。	•
Durable, lightweight aluminum shell	٠	٠
Aluminum shells covered in foam, gel and breathable, extremely stretchy material that expands with the gel to allow complete immersion	٠	٠
Geo-Matrix gel insert along the spine to reduce risk of shear and pressure sores along the spinal column	٠	٠
High-density foam core for increased support and comfort	٠	٠
Unique ability to add custom gel or foam pieces for additional customized support and contour	٠	٠
Additional layer of Geo-Matrix gel over full back for maximized comfort and immersion to further reduce risk of skin ulcers		٠
More supple immersion for kyphotic and bony prominence pressure reduction		٠
Custom height and width sizes are availablee	٠	٠

HOW TO **ORDER** YOUR **SILVERBACK** 

MAX

Μ	16	S	Α
Style	Width	Height	Contour
P-Prime	16"	S-Standard (16")	A-Active
M-Max	18"	T-Tall (18")	D-Deep
	20"	XT-Extra-Tall (20")	XD-Extra Deep
	Custom width	Custom height	

PRIME



#### Warranty Information

We proudly guarantee our Silverback line of wheelchair backs with a two-year warranty because we are committed to ensuring your complete satisfaction.

