



YOU'RE CLEARED FOR TAKE OFF

In the tradition of high-performance motors and gearboxes, the M5 takes control with a 7.5 mph motor package and superior torque. Coupled with a 14" seat elevator that maintains full suspension at any seating position and drive speed, piloting the M5 Corpus will help you take adventure to a new level.

3.2MPH @ 14"
OF ELEVATION +
FULL SUSPENSION
AT ANY SPEED,
POSITION, OR
ELEVATION







The M5 Corpus™ is unlike anything you have ever experienced. Our unique suspension geometry provides increased stability even in 14" of seat elevation. The FlexLink Pro™ suspension combines the advanced vibration suppression of FlexLink™ with independent drive wheel suspension, giving you the power to traverse all of life's obstacles while maintaining your desired seating position.

ACTIVE REACH* | REACH FOR MORE

Your world isn't one dimensional. Your seat elevator shouldn't be either. ActiveReachTM technology is the ultimate advantage in accessing your environment. By combining 14" of seat elevation and up to 20° of forward tilt, Permobil's ActiveReachTM technology enables you to reach beyond obstructions, into high cabinets, and can even assist in transfers. Truly putting your world within reach.

the . By ward up to even th.

14" OF ELEVATION + 20° FORWARD TILT = ULTIMATE ACCESSIBILITY.



COLOR & HYDROGRAPHIC OPTIONS



Radiant Red



Power Pink



Ocean Blue



Phantom Grey



Mossy Oak® Breakup® Country[™]



CDG Poly

FABRIC OPTIONS



Leatherette



Mesh

















Carbon Fiber

SPECIFICATIONS

Maximum User Weight	300 lbs
Maximum Speed	7.5 mph
Range	16 mi*
Base Width	25"
Base Length	37.5"
Minimum Turning Radius	21.5"
Ground Clearance	3"
Recommended Battery Type	Group 24
Weight Including Batteries	430 lbs
Drive Electronics	R-Net 120A
Seat-to-Floor Height	17.5", 18.5", 19.5"

Optional Seat Elevator	14"
Tilt Options	0° – 50° Posterior 5°, 10° & 20° Anterior
Recline Options	85° – 120° (manual) 85° – 180° (power)
Legrest Angle	85° – 170°
Backrest Height	20", 23" – 28" (by 1" incr.)
Seat Depth	14" – 22" (by 1" incr.)
Seat Width	17" – 23" (by 2" incr.)
Armrest Pad Lengths	10", 13", 16", 18"

[†]Specifications may vary based upon configuration. Please refer to your owners manual for warnings and instructions for the safe operation of your wheelchair in a variety of driving conditions.

^{*}Actual driving range will vary based on driving and battery conditions.